

FLORENCE: Clinician Dashboard - User Manual

1. Introduction

The FLORENCE Clinician Dashboard is built for speed and secured for healthcare. It provides medical professionals with a streamlined, high-priority command centre that turns scattered patient logs into actionable clinical insights. The goal: Less noise. More signal.

2. Dashboard & Dynamic Risk Engine (Home Screen)

2.1 Overview of the Home Screen

Upon logging in, you are presented with the main command centre. This screen is designed to highlight patients requiring immediate attention.

Clinician Dashboard

Search patient by name or ID

J

Jordan

ID: 396

Type 2 Diabetes, Stage 2 Hypertension

Last Update
May 25

B

Basil Agas

ID: 106

Type 2 Diabetes

Last Update
3h ago

Basil Agas

3h ago

High HbA1c Reading

High HbA1c: 10.0%

View Details

Jordan

May 25, 2:27 AM

High Blood Pressure

High Systolic BP: 150.0 mmHg

View Details

Jordan

May 25, 2:27 AM

High Glucose Reading

High Glucose: 10.8 mg/dL

View Details

Jordan

May 24, 2:22 AM

High Blood Pressure

High Systolic BP: 145.0 mmHg

View Details


2.2 Key Features of the Home Screen

1. **Search & Filter [Box 1]:** Quickly find patients by typing their name or ID. Tap the filter icon to expand advanced filters (e.g., filtering by specific risk levels or update recency).
 2. **Dynamic Sorting [Box 2]:** The patient list is NOT alphabetical. It is driven by a live Dynamic Risk Engine. Patients flagged as "High" or "Medium" risk automatically float to the top of your queue based on real-time continuous data flows from their patient app.
 3. **Priority Alerts [Box 3]:** If a patient breaches critical thresholds (e.g., consecutive fasting glucose readings > 10.0 mmol/L or hypertensive crisis readings), a **Priority Alert** card appears here. It displays the patient's name, the specific trigger, and the time elapsed. Tap **View Details** to investigate immediately.
-

3. Patient Overview & Demographics

3.1 Holistic Patient View

Clicking on any patient from your list opens their comprehensive profile, defaulting to the **Overview** tab.



3.2 Understanding the Overview Tab

1. **Demographics & Baseline [Box 1]:** Displays essential patient identifiers, age, gender, and emergency contact details. It also features an **Edit** button (pencil icon) to update profile details.
2. **Current Status Grid [Box 2]:** A stack of vital metrics cards showing the latest logged values. Each card is colour-coded by its current risk status:
 - **Green (Normal):** On track and within target thresholds.
 - **Amber (Elevated/Low):** Out of range; requires monitoring.
 - **Red (High/Critical):** Significantly out of range; requires immediate clinical review.

- **Grey (No Data):** No readings logged yet.

3. **Clinical Notes [Box 3]:** Displays persistent, chronological notes from you and your clinical team. Tap the **Add Note** icon to append a new entry.

4. Historical Data & Dynamic Analytics

4.1 Navigating to Historical Data

To view detailed trends and long-term patient data, navigate to the **Historical Data** tab.



Historical Data

Overview

Medical Profile

Historical Data

Body Mass Index



22.4

Normal

Weight

64.7 kg

Height

170 cm

Metric

Imperial



Glucose

10 mg/dL

Low >



HbA1c

10%


High >



Blood Pressure

80/60 mmHg

Low >



4.2 Glucose Analytics Screen

Clicking on the **Glucose** summary card in the **Historical Data** tab opens the dedicated Glucose Analytics Screen.



Glucose Analytics

1

24 Hours

✓ 7 Days

12 Months



1 Jun - 7 Jun 2026



Time in Range



0% Low

100% In Range

0% High

Average



112 mg/dL

Variability



12.5%


Readings

27

Lowest



90 mg/dL



Key Features:

- 1. Timeframe Filters [Box 1]:** Toggle between daily, weekly, and monthly views. Use the chevron arrows to navigate back and forth in time.
- 2. Smart Charts [Box 2]:** The charts dynamically scale their intervals so that even months of daily data remain readable without label overlap. Shaded bands represent the patient's personalised target safe zones.
- 3. Status Pills [Box 3]:** Below the charts, individual logged entries are listed chronologically and flagged with status pills (e.g., "ABOVE TARGET", "BELOW TARGET", "NORMAL") for rapid clinical assessment.

4.3 HbA1c Analytics Screen

Clicking on the **HbA1c** summary card in the **Historical Data** tab opens the dedicated HbA1c Analytics Screen.

24 Hours

✓ 7 Days

12 Months



1 Jun - 7 Jun 2026



Actual vs. Goal

1

10.0%



You


7.0%



Target



You are above your target of <7%



Key Features:

1. **Actual vs. Goal [Box 1]:** A clean, side-by-side bar chart comparing the patient's latest HbA1c reading against their target. A status card below immediately flags whether they are within target.
2. **HbA1c Trends [Box 2]:** A line chart showing long-term control over months. Shaded bands represent the patient's target safe zones.
3. **Status Pills [Box 3]:** Below the charts, individual logged entries are listed chronologically and flagged with status pills (e.g., "ABOVE TARGET", "BELOW TARGET", "NORMAL") for rapid clinical assessment.

4.4 Blood Pressure Analytics Screen

Clicking on the **Blood Pressure** summary card in the **Historical Data** tab opens the dedicated Blood Pressure Analytics Screen.

Overview

1

Target Ranges



Systolic	90 - 120 mmHg
Diastolic	60 - 80 mmHg

Avg Systolic ↓

108 mmHg


Avg Diastolic ↓


72 mmHg

2

Blood Pressure Trends

3






Key Features:


1. **Target Ranges [Box 1]:** Displays the patient's personalised target ranges for both Systolic and Diastolic pressure.
2. **Averages [Box 2]:** Displays the patient's average Systolic and Diastolic pressure for the selected period.
3. **Blood Pressure Trends [Box 3]:** A line chart showing both Systolic and Diastolic lines. Shaded bands represent the patient's target safe zones.
4. **Status Pills [Box 4]:** Below the charts, individual logged entries are listed chronologically and flagged with status pills (e.g., "HIGH", "NORMAL", "LOW") for rapid clinical assessment.

4.5 Cholesterol Analytics Screen

Clicking on the **Cholesterol** summary card in the **Historical Data** tab opens the dedicated Cholesterol



Analytics Screen.



Key Features:

1. **Cholesterol Ratio [Box 1]:** A clean, circular gauge comparing the patient's HDL against their Non-HDL cholesterol. A status card below immediately flags whether they are within target.
 2. **Metric Selector [Box 2]:** Toggle between Total, LDL, HDL, and Triglycerides views.
 3. **Cholesterol Breakdown [Box 3]:** A line chart showing long-term control over months. Shaded bands represent the patient's target safe zones.
 4. **Status Pills [Box 4]:** Below the charts, individual logged entries are listed chronologically and flagged with status pills (e.g., "HIGH", "NORMAL", "LOW") for rapid clinical assessment.
-

4.6 Physical Activity Analytics Screen

Clicking on the **Physical Activity** summary card in the **Historical Data** tab opens the dedicated Physical Activity Analytics Screen.

← Activity Analytics

🏃 Movement on 4 Jun 2026



30

minutes

across 1 sessions



0

Steps



180

Calories

1

📅 Activity Streak

June 2026

M

T

W

T

F

S

S

1

2

3

4

5

6

7

8

9

10

11

12

13

14

15

16

17

18

19

20

21

22

23

24

25

26

27

28


29

30

Less ■■■ More



2



Key Features:

1. **Today's Movement [Box 1]:** Displays the patient's active minutes and sessions for today.
2. **Activity Streak [Box 2]:** A 28-day heatmap grid showing a visual representation of the patient's activity. Green represents active days, light green represents less active days, and grey represents inactive days.
3. **Weekly Consistency [Box 3]:** A bar chart showing the patient's steps vs target. Green represents days where the target was met, amber represents days where the target was partially met, and red represents days where the target was not met.

4.7 Diet Log & Nutrition Summary

Clicking on the **Diet Log** summary card in the **Historical Data** tab opens the dedicated Diet Log & Nutrition Summary Screen.

 Diet Log

1



BREAKFAST

Jun 4, 12:00 PM

Oatmeal with almond milk and blueberries (1 serving)

Carbs 0g

Protein 0g

Fat 0g



LUNCH

Jun 4, 12:00 PM

Grilled chicken breast with brown rice and steamed bok choy (1 serving)

Carbs 0g

Protein 0g

Fat 0g



DINNER

Jun 4, 12:00 PM

Steamed fish with ginger and stir-fry vegetables (1 serving)

Carbs 0g

Protein 0g

Fat 0g

View all entries

 Automated Actions Log

3



Alert

Jun 8, 12:09 PM

Hypertensive Crisis: BP of 180/120 mmHg is dangerously high. Seek medical attention immediately.



Encouragement

Jun 8, 12:08 PM

Proud of your progress!: No activity logged in 2 days. Even a 10-minute walk helps manage glucose and mood!

Key Features:

1. **Diet Log List [Box 1]:** Displays the patient's logged meals chronologically. Each entry shows the meal type, timestamp, and food items.

- 2. **Nutrient Chips [Box 2]:** Displays the patient's macronutrient breakdown for each meal.
- 3. **Automated Actions Log [Box 3]:** Displays the patient's automated actions log. Each entry shows the action type, timestamp, and description.

5. Taking Clinical Action


5.1 Managing Medical Conditions & Medications

Navigate to the **Medical Profile** tab to manage the patient's clinical background, target thresholds, and active prescriptions. This tab serves as the clinical core for personalising patient care.



Detailed Button Functions & Clinical Workflows:

1. Medical Conditions Card [Box 1]




This card displays the patient's active and resolved diagnoses.

- **Add New Condition Button (Plus Icon):** Tap this button to open the **Add Medical Condition** dialog.
 - *Condition Name:* Enter the diagnosis (e.g., "Type 2 Diabetes", "Hypertension").
 - *Status Context Dropdown:* Select **Active** (for ongoing management) or **Resolved** (for historical context).
 - *Diagnosed Date:* Tap to open a calendar date picker to log the exact onset date.
 - *Add Entry Button:* Submits the record to the backend, immediately updating the patient's clinical profile.
- **Active / Resolved / All Segmented Filter:** Tap these segments to instantly filter the list. This helps you quickly isolate active issues during a consultation without being distracted by resolved history.
- **More Options Button (Three-Dot Icon):** Located next to each condition. Tap this to open a context menu:
 - *Mark as Resolved / Active:* Instantly toggles the condition's status and moves it to the appropriate filtered list.
 - *Edit Details:* Re-opens the dialog to correct names or onset dates.

- *Remove*: Permanently deletes the condition from the patient's record after a confirmation prompt.

2. Health Thresholds Card [Box 2]



The screenshot shows a card titled "Health Thresholds" with a list of vital signs and their target ranges. A red box highlights the edit icon (pencil) in the top right corner. The list includes:


Vital Sign	Target Range
Glucose	72.1 - 180.2 mg/dL
Total Cholesterol	2.0 - 10.0 mmol/L
LDL Cholesterol	0.0 - 10.0 mmol/L
HDL Cholesterol	2.0 - 10.0 mmol/L
Triglycerides	0.0 - 2.0 mmol/L
BP (Diastolic)	60.0 - 80.0 mmHg
BP (Systolic)	90.0 - 120.0 mmHg
Healthy BMI	18.5 - 24.9
HbA1c	4.0 - 7.0 %

This card displays the patient's

personalised target ranges for vitals (Glucose, Blood Pressure, HbA1c, Cholesterol, and BMI). These thresholds are critical because the **Dynamic Risk Engine** uses them to calculate the patient's risk level and float them to the top of the home screen queue.

- **Edit Button (Pencil Icon)**: Tap this to open the **Health Thresholds** management dialog.
 - *Threshold List*: Displays all vitals with their current min/max values.
 - *Edit Single Threshold (Pencil Icon)*: Tap the edit icon next to any vital (e.g., Glucose) to open the **Edit Threshold** dialog.
 - *Minimum / Maximum Value Fields*: Enter the new clinical targets. The units automatically match your preferred clinician settings (e.g., mmol/L or mg/dL).
 - *Save Button*: Securely updates the targets via the REST API. The patient's app and the risk engine sync instantly.

3. Current Medications Card [Box 3]





This card displays the patient's active prescriptions and schedule.

- **Add New Medication Button (Plus Icon):** Tap this to open the comprehensive **Add Medication** form.
 - *Medication Name (Smart Autocomplete):* Start typing a brand or generic name. The field queries the global medication dictionary in real-time to prevent spelling errors and ensure clinical accuracy. You can also type a custom name if it is not in the dictionary.
 - *Amount & Type:* Enter the dosage quantity (e.g., "1", "500") and select the form (e.g., "Tablet", "Capsule", "Injection", "ml").
 - *Frequency Dropdown:* Select how often the medication should be taken (e.g., "Twice daily", "Three times daily").
 - *Specific Timings (Dynamic Fields):* The form dynamically generates dropdown fields based on your selected frequency (e.g., showing "1st Dose" and "2nd Dose" for twice-daily frequency). Select precise clinical instructions for each dose (e.g., "BEFORE_BREAKFAST", "AFTER_DINNER", "BEFORE_BED").
 - *Add Medication Button:* Saves the prescription. This automatically generates a daily logging schedule in the patient's mobile app.
- **More Options Button (Three-Dot Icon):** Located next to each medication. Tap this to:
 - *Edit Details:* Modify dosage, frequency, or timing instructions.
 - *Remove:* Discontinue and delete the medication from the patient's active schedule.

5.2 Adding Clinical Notes

To record observations, treatment adjustments, or consultation summaries, tap the **Add Note** button on the patient's Overview tab.





Clinical Actions:

1. **Enter Note [Box 1]:** Type your clinical observations or treatment adjustments.
2. **Save Note [Box 2]:** Tap **Save Note**. Because the dashboard operates on a strict "Zero Direct Database" policy, your note is securely routed through the FLORENCE middleware REST API. The note immediately syncs across the platform, updating the priority queue and context for your entire clinical team in real-time.

6. Clinician Profile & Preferences

Accessing the Profile

To view your account details, update your personal information, or modify your unit preferences, tap the **Profile** icon in the top-right corner of the Home Screen.




Image Highlight Instructions (No image file, text-only directions for UI layout):

Please capture a screenshot of the **Clinician Profile** screen and overlay **four red highlight boxes**, numbered as follows:

- **[Box 1]** Around the **Account Information Card** (containing your name, gender, and mobile number) with the **Edit Profile** button.
- **[Box 2]** Around the **Unit Preferences Card** (containing Glucose Unit and Cholesterol Unit selectors).
- **[Box 3]** Around the **System Information Card** (containing your email) and the **Organisation Information Card** (containing uneditable organisation details: name, email, and phone number).
- **[Box 4]** Around the **Logout Button** at the bottom.

Detailed Card Descriptions & Clinical Actions:

1. Account Information Card [Box 1]

←

Clinician Profile

DS

Dr. Sarah

Organisation ID: 4

Account Information

Full Name

Dr. Sarah

Gender

Female

Mobile Number


0138200334

Edit Profile

This card displays your personal clinical identity details used across the platform.

- **What it is:** A secure card containing your editable profile details: **Full Name**, **Gender**, and **Mobile Number**.
- **What it does:** Allows you to keep your contact details and identity up to date so that patients and other clinical staff can identify and contact you correctly.
- **How to use it:**
 1. Tap the **Edit Profile** button in the bottom-right corner of the card. The fields will switch from read-only to editable text fields.
 2. Modify your **Full Name** or **Mobile Number** as needed.
 3. Select your **Gender** from the dropdown menu.
 4. Tap **Save Changes** to commit the updates to the backend database via the REST API, or tap **Cancel** to discard your edits and restore the previous values.

2. Unit Preferences Card [Box 2]



This card controls how clinical metrics are displayed to you throughout the entire dashboard.

- **What it is:** A preference card containing selectors for **Glucose Unit** and **Cholesterol Unit**.
- **What it does:** Customises the display units for all patient vitals across the dashboard. Changing these units automatically converts and displays all historical readings, charts, and thresholds in your preferred format (e.g., converting glucose from mmol/L to mg/dL) without altering the underlying base data in the database.
- **How to use it:**
 1. Tap the **Glucose Unit** row to open a bottom sheet selector. Choose between **mmol/L** and **mg/dL**.
 2. Tap the **Cholesterol Unit** row to open a bottom sheet selector. Choose between **mmol/L** and **mg/dL**.
 3. The platform instantly recalculates and updates all charts, tables, and status cards across the app. The preferences are saved to your profile via the REST API and persist across logins.

3. System Information & Organisation Information Cards [Box 3]

System Information

Email Address



clinician@test.user

Organisation Information

Organisation Name



Metropolis Health

Organisation Email



admin@metropolicehealth.com

Organisation Phone Number



+6 082-507 333

These cards display your

immutable system credentials and detailed organisational context.

- **What they are:** Two read-only cards. The **System Information Card** displays your registered **Email Address**. The **Organisation Information Card** displays your assigned facility's **Organisation Name**, **Organisation Email**, and **Organisation Phone Number**.
- **What they do:** Provide essential system and facility context. Your email is your unique login identifier, while the organisation details link your account to your specific healthcare facility, determining which patients, priority alerts, and facility contact channels you have permission to access.
- **How to use them:** These sections are strictly read-only for security, compliance, and audit purposes. If your email or organisation assignment needs to be changed, please contact your system administrator.

4. Logout Button

This button allows you to securely terminate your active session.

- **What it is:** A red-outlined button at the bottom of the profile screen.
- **What it does:** Instantly invalidates your active Supabase authentication session, clears cached patient data from the device memory to prevent unauthorised access, and returns you to the secure login screen.
- **How to use it:**
 1. Tap the **Logout** button.
 2. A confirmation dialog will appear asking: "Are you sure you want to sign out?"
 3. Tap **Sign Out** to confirm. You will be immediately redirected to the login screen.